DCI™ Dynamic Cervical Implant

Package contents

The package contains one Dynamic Cervical Implant (DCI™).

Description

The DCI™ implant stabilizes the cervical spine providing stable, controlled motion. Flexion and extension are allowed relative to a neutral position.

Material

The DCI™ implant is made from wrought titanium 6-aluminum 4-vanadium alloy (ISO 5832-3). It is supplied sterile and is available in numerous sizes.

Intended use

The DCI™ implant is intended for permanent implantation between the anterior cervical disc space at 1 to 3 levels from C3 to C7 in skeletally mature patients and controls segmental motion in cases of cervical disc herniation, cervical degenerative disc disease (DDD) and cervical canal stenosis (central or foraminal) with or without myeloradiculopathy.

Indications

The DCI™ implant is indicated for anterior implantation into the cervical disc space at 1 to 3 levels from C3 to C7 and controls segmental motion in cases of cervical disc herniation, cervical DDD and cervical canal stenosis (central or foraminal) with or without myeloradiculopathy in patients with or without neck pain.

Contraindications

Do not use the DCI™ implant in cases of:

- Any medical or surgical condition precluding the potential benefit of spinal surgery
- Acute or chronic systemic, spinal or localized infections
- Uncontrolled systemic and metabolic diseases
- Obesity
- Pregnancy
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Lack of patient cooperation
- Foreign body sensitivity to the implant material
- Severe mechanical instability
- Osteoporosis
- Significant osteopenia
- Vertebral fractures and/or tumors and at any cervical level that has absence of motion on preoperative flexion/extension radiographs

Adverse reactions

Adverse reactions may include:

- Clinical failure (i.e. pain or injury) due to bending, loosening, wear and tear, fracture of implant, loss of fixation, dislocation and/or migration
- Pain and/or abnormal sensations due to the presence of the implant
- Primary and/or secondary infections
- Allergic reactions to implant material
- Neurological injury
- Vertebrae fracture
- Injury to vessels, nerves and organs
- Hematoma and/or impaired wound healing; hemorrhage
- Venous thrombosis, lung embolism, and cardiac arrest
- Death

Safety precautions

- Prior to use, thoroughly read these instructions for use and become familiar with the surgical technique.

- Keep the instructions for use accessible to all staff.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques. Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrect operating techniques, the limitations of treatment methods or inadequate asepsis.
- Under no circumstances may modular implant components from different suppliers be combined.
- Each patient's record shall document the implant used (name, article number, lot number).
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed about post-surgical behavioral requirements.
- Damage to the weight-bearing structures can give rise to loosening, dislocation and migration, as well as other complications. To ensure the earliest possible detection of such catalysts of implant dysfunction, the implant must be checked periodically post operative using appropriate techniques.

- Never reuse an implant. Although the implant may appear undamaged, previous stresses may have created non-visible damage that could result in implant failure.

- Never use implants if the packaging is damaged. An implant with damaged packaging might be damaged itself and thus may not be used.

- Never use implants that are past their expiration date.

MRI compatibility

Non-clinical testing has demonstrated that the DCI™ Implant is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5-Tesla (1.5T) or 3.0-Tesla (3.0T)
- Spatial gradient field of up to:
 - 9,760 G/cm (97.60 T/m) for 1.5T systems
 - 4,880 G/cm (48.80 T/m) for 3.0T systems.
- Maximum whole body averaged specific absorption rate (SAR) of up to:
 - 4.0 W/kg for 15 minutes of scanning in Normal Operating/First Level Modes at 1.5T
 - 4.0 W/kg for 15 minutes of scanning in Normal Operating/First Level Modes at 3.0T.

3.0T RF heating

In non-clinical testing with body coil excitation, the DCI™ Implant produced a temperature rise of less than 1.5 °C at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 3.0T Siemens Trio (MRC20587) MR scanner with SYNGO MR A30 4VA30A software.

1.5T RF heating

In non-clinical testing with body coil excitation, the DCI™ Implant produced a temperature rise of less than 2.0 °C at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 1.5T Siemens Espree (MRC30732) MR scanner with SYNGO MR B17 software.

The RF heating behavior does not scale with static field strength. Devices which do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

MR Artifact

In testing using a 3.0T system with spin-echo sequencing, the shape of the image artifact follows the approximate contour of the device and extends radially up to 23 mm from the implant.

Storage, inspection and sterilization

Storage

The implant is individually packed in protective packaging that is labeled according to its contents.

 The implant is sterilized with gamma sterilization (25 kGy minimum).

- Always store the implant in the original protective packaging.
- Do not remove the implant from the packaging until immediately before use.
- Store the implant in a dry and dust-free place (standard hospital environment).

Disinfection / cleaning

The DCI™ implant is not designed to be disinfected or cleaned by the user.

Restoration

The DCI™ implant is not designed to be re-sterilized by the user.

Procedure

The DCI™ implant must be implanted only with the applicable DCI™ instruments. The DCI™ instruments are available from the manufacturer at any time.

Preoperative

The operating surgeon draws up an operation plan specifying and documenting the following:

- Implant component(s) and their dimensions.
- Proper position of the implant component(s) between the anterior cervical disc space.
- Determination of intraoperative orientation points.

The following conditions must be fulfilled prior to application:

- All required implant component(s) are readily available.
- Highly aseptic operating conditions are present.
- All requisite implantation instruments must be available and in working order, including the DCI™ trial

Never use or process damaged or defective instruments. Contact your local representative or dealer for repair or replacement.

The use of an instrument for tasks other than those for which they are intended may result in damaged / broken instruments or patient injury.

- The operating surgeon and operating room team must be thoroughly familiar with the operating technique, as well as the range of implants and instruments to be applied. Complete information on these subjects must be readily available at the workplace.
- The operating surgeon must be especially trained in spinal surgery, biomechanical principles of the spine and the relevant operating techniques.

The operating procedure must be explained to the patient, and the patient's understanding of the following information must be documented:

- The patient is aware of the risks associated with neurosurgery, general surgery, orthopedic surgery and with general anesthesia.

- The patient has been informed about the advantages and disadvantages of the implant procedure and about possible alternative treatments.
- The implant can fail due to excessive load, wear and tear, or infection.
- The service life of the implant is determined by body weight and physical activity. The implant must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may be necessary if the implant fails.
- The patient must have their physician carry out followup examinations of the implant at regular intervals.

Intraoperative

Prior to use, please read the DCI™ product brochure and become familiar with the surgical technique, the corresponding implants and the instruments.

- Prior to use, verify the integrity of the sterile packaging. Never use implants if the packaging is damaged.

- Prior to use, check the product expiration date. Never use implants that are past their expiration date.

The disc is excised and the intervertebral space is prepared for implantation. Microsurgical decompression is performed when necessary. The intervertebral space is measured utilizing the DCI™ trial with sleeve and turning knob. Load the DCI™ implant on the insertion instrument or inserter for insertion and final positioning.

Postoperative

- Reiterate preoperative instructions to the patient.
- Ensure that the patient is aware of physical activity restrictions and possible adverse reactions.

Revision surgery / implant removal

The DCI™ implant is intended for permanent implantation and is not intended for removal. However, removal may be advisable in the following situations:

- Implant breakage
- Pain due to the implant
- Infection

Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use. Technical alterations reserved.




For further information

Please contact Paradigm Spine or your authorized representative if further information on this product is needed.

Patent

For patent information see www.paradigmspine.com

Symbols

	Manufacturer
	Date of manufacture
	Use-by date
	Lot number
	Catalog number
	Content of usable units(s)
	Sterilized using irradiation
	Do not resterilize
	Do not use if package is damaged
	Keep dry
	Do not reuse
	Consult instructions for use at this website.
	Caution
	MR conditional
	CE-marking according to directive 93/42/EEC

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