## Streamline® TL/MIS Fenestrated Spinal Fixation System

### Package contents

The package contains one or more of the various components of the Streamline TL/MIS fenestrated Spinal Fixation System:

- Streamline TL/MIS fenestrated screw Ø 5,5
- Streamline TL/MIS fenestrated screw Ø 6,5
- Streamline TL/MIS fenestrated screw Ø 7,5
- Streamline TL/MIS fenestrated screw Ø 8,5
- Streamline TL/MIS fenestrated screw Ø 9,5
- Streamline TL/MIS screw Ø 4,5 non fenestrated

The content of each package is evident from the respective product label.

### Description

The modular Streamline TL/MIS Spinal Fixation System is used for single or multi-segmental posterior rigid stabilization of the thoracolumbar spine. The system is based on fixation with polyaxial pedicle screw, fenestrated and allows the application with bone cement.

Streamline TL/MIS screw Ø 4,5 non fenestrated is not intended for the application with bone cement.

The Streamline TL/MIS, fenestrated screw is intended for the use with HPS™ rods and RTI Surgical Streamline TL/MIS rods.

If the Streamline TL/MIS, fenestrated screw is used in combination with RTI Streamline TL/MIS Implants and RTI Streamline instruments this instruction for use (IFU) is primarily applicable.

### Material

All Streamline TL/MIS pedicle screws, fenestrated are made from titanium 6-aluminum 4-vanadium alloy (ISO 5832-3).

The used material contains no nickel.

Not made with natural rubber latex.

### Intended use

The Streamline TL/MIS Spinal Fixation System is a pedicle screw-based implant system for posterior rigid stabilization of the thoracolumbar spine in skeletally mature patients. The Streamline TL/MIS is intended for long-term implantation and for single-use only. The Streamline TL/MIS pedicle screws, fenestrated are delivered sterile and are not intended for re-processing by the user or third parties. The Streamline TL/MIS Spinal Fixation System shall only be used by surgeons who are trained and familiar with the implant components, instruments, and surgical technique.

### Indications

#### Posterior Rigid Stabilization

The Streamline TL/MIS Fenestrated Spinal Fixation System is intended for posterior rigid stabilization of the thoracolumbar spine as adjunct to fusion for the following indications:

- degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies);
- spondylolisthesis;
- trauma (i.e., fracture or dislocation);
- spinal stenosis;
- curvatures (i.e., kyphosis, and/or lordosis);
- tumor;
- pseudoarthrosis;
- and failed previous fusion

### Application with bone cement

The Streamline TL/MIS fenestrated screw, is indicated for the use with bone cement in patients with diminished bone quality (e.g. osteoporosis, osteopenia, metastatic disease) and shall only be used with components for posterior rigid stabilization of the Streamline TL/MIS Spinal Fixation System for mono- or multisegmental stabilization of the thoracic and lumbar spine for supporting fusion.

Streamline TL/MIS screw Ø 4,5 non fenestrated is not intended for the application with bone cement.

### Contraindications

#### General contraindications

Generally do not use the Streamline TL/MIS Spinal Fixation System for the following conditions:

- All medical and surgical conditions that exclude the advantages of spine surgery
- Damage to the bone structure that prevents a stable implantation
- Acute or chronic systemic, spinal or localized infections
- Impaired wound healing
- Active systemic and metabolic diseases
- Obesity (BMI ≥30)
- Fever
- Pregnancy
- Medication and substance abuse, alcoholism
- Lack of patient cooperation
- Mental illness
- Sensitivity to foreign bodies such as implant materials
- Serious osteopenia if used without bone cement
- Osteoporosis if used without bone cement
- Neuromuscular disorders or illnesses
- Bone tumors in the proximity of the implant fixation
- Fracture
- Other cases that are not listed among the "Indications"

#### Contraindications with use of bone cement

- Poor visualizing at Fluoroscopy
- Thrombophilia
- Patients with severe heart and/or lung failure
- Patients with known sensitiveness to elements of the bone cement

#### Refer to the IFU of the bone cement supplier for further contraindications

### Adverse reactions

Adverse reactions may include:

- Clinical failure (i.e. pain or injury) due to bending, loosening, wear and tear, fracture of implant, loss of fixation, dislocation and/or migration of the implant
- Pain and/or abnormal sensations due to the presence of the implant
- Primary or secondary infections
- Allergic reactions to implant material
- Neurological injury
- Vertebrae fracture
- Injury to vessels, nerves, the spinal cord, and organs
- Hematoma and/or impaired wound healing; hemorrhage
- Venous thrombosis, lung embolism, and cardiac arrest
- Lack of fixation

- No or delayed fusion
- Change to the spinal column curvature
- Impairment of the gastrointestinal, urological and/or reproductive system
- Pain or malaise
- Bursitis
- Decrease of bone density due to avoiding load
- Osteoporosis or fracture above or below the affected area of the spine
- Persistence of the symptoms to be treated with the implantation
- Death

#### Adverse reaction with use of bone cement

- Hypersensitivity reactions by susceptible patients which may cause to an anaphylactic shock
- Tissue damage, nerve damage or circulatory problems by cement leakage
- Micro movement of the cement against the surface the bone due to inadequate fixation
- Fat embolism
- decrease in blood pressure
- myocardial infarction
- cardiac arrest
- cerebrovascular accident
- pulmonary embolism
- cardiac embolism
- short-term cardiac conduction disorder
- sudden death

#### Safety precautions

- Prior to use, thoroughly read these instructions for use and become familiar with the surgical technique.
- If the Streamline TL/MIS pedicle screw, fenestrated is used in combination with RTI Streamline TL/MIS Implants and RTI Streamline instruments this IFU is primarily applicable

- Keep the instructions for use accessible to all staff.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrect operating techniques, the limitations of treatment methods or inadequate asepsis.
- Under no circumstances may modular implant components from different suppliers be combined.
- Damage to the weight-bearing structures can give rise to loosening, dislocation and migration, as well as other complications.
- To ensure the earliest possible detection of such catalysts of implant malfunction, the implant must be checked periodically post-operatively using appropriate techniques.

- Never reuse an implant. Although the implant may appear undamaged, previous stresses may have created non-visible damage that could result in implant failure.

- Never use implants if the packaging is damaged. An implant with damaged packaging might be damaged itself and thus may not be used.

- Never use implants that are past their expiration date.

#### MRI compatibility

The Streamline TL/MIS spinal fixation system has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Streamline TL/MIS spinal fixation system in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

#### Storage, inspection and sterilization

##### Storage

The implant is individually packed in protective packaging that is labeled according to its contents.

**STERILE R** The implant is sterilized with gamma sterilization (25 kGy minimum).

- Always store the implant in the original protective packaging.
- Do not remove the implant from the packaging until immediately before use.

- Store the implant in a dry and dust-free place (standard hospital environment).

#### Disinfection / cleaning

The Streamline TL/MIS pedicle screws fenestrated are not designed to be disinfected or cleaned by the user.

#### Resterilization

- The Streamline TL/MIS pedicle screws fenestrated are not designed to be re-sterilized by the user.

#### Technique

The Streamline TL/MIS Fenestrated Spinal Fixation System must be implanted only with the applicable RTI surgical TL/MIS instruments and Paradigm Spine Streamline TL/MIS instruments. The instruments are available from the manufacturer at any time.

#### Preoperative

The operating surgeon draws up an operation plan specifying and documenting the following:

- Implant component(s) and their dimensions.
- Proper positioning of the implant components.
- Determination of intra-operative orientation points. The following conditions must be fulfilled prior to application:
- All required implant component(s) are readily available.
- Highly aseptic operating conditions are present.
- All requisite implantation instruments must be available and in working order.

*Never use or process damaged or defective instruments. Contact your local representative or dealer for repair or replacement.*

*The use of an instrument for tasks other than those for which they are intended may result in damaged / broken instruments or patient injury.*

- The operating surgeon must be especially trained in spinal surgery, biomechanical principles of the spine and the relevant operating techniques.

- The operating surgeon and operating room team must be thoroughly familiar with the operating technique, as well as the range of implants and instruments to be applied. Complete information on these subjects must be readily available at the workplace.

The operating procedure must be explained to the patient, and the patient's understanding of the following information must be documented:

- The patient is aware of the risks associated with neurosurgery, general surgery, orthopedic surgery and with general anesthesia.
- The patient has been informed about the advantages and disadvantages of the implant procedure and about possible alternative treatments.
- The implant can fail due to excessive load, wear and tear, or infection
- The service life of the implant is determined by body weight and physical activity. The implant must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may be necessary if the implant fails. The patient must have their physician carry out follow-up examinations of the implant at regular intervals.

#### Intra-operative

Prior to use, please read and become familiar with the Streamline TL/MIS Spinal Fixation System surgical technique, the corresponding implants and the instruments.

- Prior to use, verify the integrity of the sterile packaging. Never use implants if the packaging is damaged.

- Prior to use, check the product expiration date. Never use implants that are past their expiration date.

#### Posterior Rigid Stabilization

- Only bend rods with the intended instrument.
- Never bend back pre-bent rods.

#### Application with bone cement

If the use of bone cement is required, it has to be introduced via the Streamline TL/MIS Polyaxial pedicle screw, fenestrated utilizing the Streamline TL/MIS Cement injection set. The bone cement must have a medium viscosity and is intended only for posterior rigid stabilization.

Streamline TL/MIS screw Ø 4,5 non fenestrated is not intended for the application with bone cement.

#### General information regarding technique

- Do not use implants that are damaged with scratches, notches or unintentional bending.
- Only use the instruments intended for the respective step.
- Only use the intended instruments (e.g., rod bender) to bend the rods.
- Use a torque limiting-handle with a counter torque instrument to avoid bending and movement of the pedicle screws during tightening.
- The proper tightening torque is reached when until the torque-limiting handle emits an audible and/or tangible "click". Apply no more torque.
- Use a distraction and/or compression instrument to obtain sagittal alignment of the implant.

#### NOTE

Detailed information regarding the use of the Streamline TL/MIS Spinal Fixation System can be found in the Streamline TL/MIS surgical technique. Further information regarding the procedure is available in workshops, product training or individual consultation with RTI Surgical or Paradigm Spine.

#### Postoperative

- Reiterate preoperative instructions to the patient.
- Ensure that the patient is aware of physical activity restrictions and possible adverse reactions.

#### Revision surgery / implant removal

The Streamline TL/MIS spinal fixation system is intended for permanent implantation and is typically not removed. However, removal of the implant can be necessary in the following situations:

- Implant breakage
- Pain due to the implant
- Infection
- Pseudarthrosis
- Allergic reactions

The surgeon shall determine whether new pedicle screws are used or if the fusion component should be placed over the existing screws. If it is necessary to replace the rod and set screws it is advisable to replace them on both sides of the spinal segment. Implants which appear to be intact may have invisible damage.

#### Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with these instructions for use. Subject to technical modifications.

#### For further information


Please contact RTI Surgical, Paradigm Spine or your medical product representative if further information is needed.

#### IMPORTANT NOTE

These instructions for use are complete at the time of publishing. The Streamline TL/MIS Spinal Fixation System may be supplemented with new products in the future. If a new implant is used which is not mentioned in this document, please read the respective instructions for use.

#### Symbols

	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Content of usable units(s)

	Sterilized using irradiation
	Do not resterilize
	Do not use if package is damaged
	Keep dry
	Do not reuse
	Consult instructions for use at this website.
	Caution
	CE-marking according to directive 93/42/EEC

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