K-Wires

Package contents

The package contains one of the following single-use instruments:

- HPS 2.0 K-wire

The content of each package is evident from the respective product label.

Description

The HPS™ 2.0 K wires are designed for use in posterior spinal stabilization procedures to guide the pedicle screw placement into the vertebral body.

Material

HPS 2.0 K-Wire are made from Nitinol (Nickel-Titanium Shape Memory Alloy).

Intended use

The HPS™ 2.0 K wires may only be used according to the directions described herein, and the instructions for use enclosed with surgical systems.

Prior to use, make sure that the instruments are sterile and have not been damaged. Do not reuse HPS™ 2.0 K wires.

Surgical Technique for posterior stabilization with the MIS approach: A longitudinal incision about 1.5 cm is made through the skin and fascia. Insert the Jamshidi Needle through incision and dock the tip on the bony anatomy of the desired level. Confirm the position by using lateral fluoroscopy. Using A-P fluoroscopy, advance the Jamshidi Needle to the pedicle entry point at the intersection of the facet and transverse process. Confirm that the tip of the Jamshidi Needle is at the center of the lateral border of the pedicle on an A-P image. Gently tap the Jamshidi Needle to engage the trocar tip into the pedicle.

Advance the Jamshidi needle through pedicle using A-P fluoroscopy to direct the tip towards the center of the pedicle. The needle should not advance more than three quarters of the way across the pedicle, starting from the lateral edge of the pedicle. Continue advancement until the needle enters the vertebral body. Confirm the placement with A/P and lateral fluoroscopy to ensure that the Jamshidi Needle does not breach the wall of the pedicle. Remove the inner stylet of the Jamshidi Needle. Select an HPS 2.0 K-Wire. Insert the HPS 2.0 K-Wire into the Jamshidi Needle. Advance the HPS 2.0 K-Wire beyond the tip of the Jamshidi Needle (approximately 20mm) to ensure an adequate fixation into the cancellous bone. Confirm the placement with A/P and lateral fluoroscopy to ensure that the HPS 2.0 K-Wire does not breach the wall of the pedicle or the vertebral body. Use caution when placing the guidewire. Markings on the wire designate 5mm increments and can be used to determine penetration depth. Additionally, the depth markers can be used to monitor unintentional guidewire advancement or rotation. Once the guidewire is placed to the desired depth, carefully remove the Jamshidi Needle while holding the guidewire.

Indications

Spinal Posterior Stabilization

The HPS™ 2.0 K wires are designed for use in spinal stabilization procedures to guide the pedicle screw placement into the vertebral body.

Contraindications

Since the HPS 2.0 K-wire are used as accessory for the applications mentioned before, there are no contraindications specific to these products.

Adverse reactions

All known reactions to orthopedic treatment have to be considered carefully. Adverse events may include but are not limited to:

- Bone resorption and poor bone formation
- Device failure
- Infection
- Tissue damage as a result of surgical trauma

Safety precautions

- The HPS™ 2.0 K wires are supplied in sterile packaging. The sterility information is indicated on the device labels.
- Do not use the product if the sterile package is damaged.
- Improper use of the HPS™ 2.0 K wires can lead to damage to the tissue, destruction of the device components, and injury to the operator, patient, or third parties.
- Avoid excessive contact pressure. This can result in damage to the device. At the same time, extreme heat is generated, which increases the risk of thermal necrosis. Excessive contact pressure can also result in breakage. Avoid jamming of the device during operation. Do not use it as a lever
- Single-use devices must be disposed of in compliance with all applicable local, state, and federal laws and regulations concerning medical waste.

- Never reuse an HPS™ 2.0 K wires
- Never use HPS™ 2.0 K wires if the packaging is damaged. An HPS™ 2.0 K wire with damaged packaging might be damaged itself and thus may not be used.
- Never use HPS™ 2.0 K wires that are past their expiration date.

Storage, inspection and sterilization

Storage

The HPS™ 2.0 K wire is individually packed in protective packaging that is labeled according to its contents.

STERILE The HPS™ 2.0 K wire is sterilized with gamma sterilization (25 kGy minimum).

- Always store the HPS™ 2.0 K wire in the original protective packaging.
- Do not remove the HPS™ 2.0 K wire from the packaging until immediately before use.
- Store the HPS™ 2.0 K wire in a dry and dust-free place (standard hospital environment).

Resterilization

The HPS™ 2.0 K wires are not designed to be re-sterilized by the user.

Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with these instructions for use. Subject to technical modifications.

For further information

Please contact RTI Surgical, Paradigm Spine or your medical product representative if further information is needed.

Symbols

	Manufacturer
	Date of manufacture
	Use-by date
LOT	Batch code
REF	Catalogue number
QTY	Content of usable units(s)
STERILE R	Sterilized using irradiation
	Do not resterilize
	Do not use if package is damaged
	Keep dry
	Do not reuse
	Consult instructions for use at this website.
	Caution
CE 0123	CE-marking according to directive 93/42/EEC

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