Streamline® TL/MIS Cement injection set

Package contents

The package contains a Streamline TL/MIS Cement injection set for the Streamline TL/MIS Fenestrated Spinal Fixation System. The Streamline TL/MIS Cement injection set consists of the following components:

- 2 x Streamline TL/MIS Cement injection kits. Each kit consists of:
 - 1 x Cement injection tubes
 - 1 x Cement injection plungers
 - 1 x Adapters for Polyaxial pedicle screw, TL
 - 1 x Adapters for Polyaxial pedicle screw, TL reduction and MIS
- 1 x Counter torque

Description

The Streamline TL/MIS Cement injection set is a sterile, single use medical aid for applying bone cement to vertebral bodies using the cement adapter and the Streamline TL/MIS Fenestrated Spinal Fixation System.

The Streamline TL/MIS Cement injection kit has a capacity of 1.5 ml.

Material

The Streamline TL/MIS Cement injection set is made from X5CrNi 18-10 austenitic Stainless Steel (1.4301) and X8CrNiS 18-9 austenitic Stainless Steel (1.4305).

Indications for use

The Streamline TL/MIS Cement injection set is intended to be used with the Streamline TL/MIS Fenestrated Spinal Fixation System, and bone cement in order to apply the bone cement into the Streamline TL/MIS Fenestrated screw, to improve the fixation of the screw within the bone. For this purpose only bone cement with a medium viscosity shall be used. Refer to the applicable IFU of the Streamline TL/MIS Fenestrated Spinal Fixation System and bone cement for detailed information and instructions for use.

Contraindications

Contraindications include the contraindications of the Streamline TL/MIS Fenestrated Spinal Fixation System and the contraindications of the used bone cement. Refer to the applicable IFU of the Streamline TL/MIS Fenestrated Spinal Fixation System and the bone cement supplier.

Adverse reactions

Adverse reactions may include:

- Clinical failure (i.e. pain or injury) due to fracture of instrument
- Primary and/or secondary infections
- Allergic reactions to instrument material
- Allergic reactions to bone cement material
- Injury to vessels, nerves, the spinal cord, and organs
- Hematoma and/or impaired wound healing; hemorrhage
- Pain or malaise
- Bursitis
- Cement leakage into the adjacent tissue

Safety precautions

- Prior to use, thoroughly read these instructions for use and become familiar with the surgical technique.

- Keep the instructions for use accessible to all staff.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.

- Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrect operating techniques, the limitations of treatment methods or inadequate asepsis.
- Under no circumstances may modular implant components from different suppliers be combined.
- Damage to the weight-bearing structures can give rise to loosening, dislocation and migration, as well as other complications.
- To ensure the earliest possible detection of such catalysts of implant malfunction, the implant must be checked periodically post-operatively using appropriate techniques.

- The Streamline TL/MIS Cement injection set is designed for single use on one patient. It may not be reused, re-sterilized or re-processed. Reusing, re-sterilizing or reprocessing may lead to the structural integrity of the product being compromised or to the product malfunctioning, which could lead to injury or illness of the patient. Furthermore, reusing, re-sterilizing or reprocessing may lead to an increased risk of contamination of the product and/or to the patient becoming infected, including but not limited to passing on infectious diseases from one patient to another.

- Never use the Streamline TL/MIS Cement injection set if the packaging is damaged. The Streamline TL/MIS Cement injection set with damaged packaging might be damaged itself and thus may no longer be used.

- Never use a Streamline TL/MIS Cement injection set that is past its expiration date.

- Each patient's record shall document the Streamline TL/MIS Cement injection set used (name, article number, lot number).
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed about post-surgical behavioral requirements.
- Caution: The operating manual of the cement must also be observed in every situation.
- The amount of cement used must be chosen by the surgeon depending on the individual anatomy of the patient and the cement used.

Storage, inspection and sterilization

Storage

The Streamline TL/MIS Cement injection set is packed in protective packaging that is labeled according to its contents.

STERILE R The instrument is sterilized with gamma sterilization (25 kGy minimum).

- Always store the instrument in the original protective packaging.
- Do not remove the instrument from the packaging until immediately before use.
- Store the instrument in a dry and dust-free place (standard hospital environment).

Disinfection / cleaning

The Streamline TL/MIS Cement injection set is not designed to be disinfected or cleaned by the user.

Resterilization

The Streamline TL/MIS Cement injection is not designed to be re-sterilized by the user.

Technique

The Streamline TL/MIS Cement injection set may only be used with the Streamline TL/MIS Fenestrated Spinal Fixation System.

1. The Adapter for Polyaxial pedicle screw, TL is screwed onto the already inserted Polyaxial pedicle screw, TL or the Adapter for Polyaxial pedicle screw, TL reduction and MIS is screwed onto the already inserted Polyaxial pedicle screw, TL reduction and MIS, by means of the Streamline TL/MIS adapter inserter and Counter torque.
2. The Cement injection tube gets filled with the prepared bone cement by a standard disposable syringe. The bone cement is applied via Cement injection plunger.
3. Hardening of the bone cement.
4. Removing the cement injection tube and plunger
5. Removing the Adapter with the help of the Counter torque and the Streamline TL/MIS adapter inserter.

Detailed information regarding the use of the Streamline TL/MIS Cement injection set can be found in the Streamline TL/MIS Fenestrated Spinal Fixation System surgical technique.

Preoperative

Refer to the applicable IFU of the Streamline TL/MIS Fenestrated Spinal Fixation System and bone cement for detailed information. for detailed information.

Never use or process damaged or defective instruments. Contact your local representative or dealer for repair or replacement.

The use of an instrument for tasks other than those for which they are intended may result in damaged / broken instruments or patient injury.

- The operating surgeon must be especially trained in spinal surgery, biomechanical principles of the spine and the relevant operating techniques.
- The operating surgeon and operating room team must be thoroughly familiar with the operating technique, as well as the range

Intra-operative

Prior to use, please read the Streamline TL/MIS and Streamline TL/MIS Fenestrated Spinal Fixation System surgical technique and become familiar with it the corresponding implants and the instruments.

- Prior to use, verify the integrity of the sterile packaging. Never use instruments if the packaging is damaged.

- Prior to use, check the product expiration date. Never use instruments that are past their expiration date.

General information regarding technique

- Do not use instruments that are damaged with scratches, notches or unintentional bending.
- Only use the instruments intended for the respective step.

NOTE

Detailed information regarding the use of the Streamline TL/MIS Cement injection set can be found in the Streamline TL/MIS Fenestrated Spinal Fixation System surgical technique. Further information regarding the procedure is available in workshops,

product training or individual consultation with RTI Surgical or Paradigm Spine.

Postoperative

- Refer to the applicable IFU of the Streamline TL/MIS Fenestrated Spinal Fixation System and bone cement for detailed information.

Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with these instructions for use. Subject to technical modifications.

For further information

Please contact RTI Surgical, Paradigm Spine or your medical product representative if further information is needed.

IMPORTANT NOTE

These instructions for use are complete at the time of publishing. The Streamline TL/MIS Fenestrated Spinal Fixation System may be supplemented with new products in the future. If a new instrument is used which is not mentioned in this document, please read the respective instructions for use.

Symbols

	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Content of usable units(s)
	Sterilized using irradiation
	Do not resterilize
	Do not use if package is damaged
	Keep dry
	Do not reuse
	Consult instructions for use at this website.
	Caution
	CE-marking according to directive 93/42/EEC

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